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STAR_android

New patient

Slide your finger up/down to scroll the screen up/down.

Patient name

Surname

First name

Hospital ID

Age years

Gender ☒ Male ☐ Female

Body frame size ☐ Small ☒ Medium ☐ Large

Blood glucose

Current mmol/L

Tab Q W E R T Y U I O P
1234 5 6 7 8 9 0
A S D F G H J K L
Z X C V B N M , . /
17:35

STAR_android

Current mmol/L

BG target range

☒ 4.4 - 8.0 mmol/L Recommended for most patients

☐ 6.0 - 9.0 mmol/L May be recommended for known Type II diabetes or other conditions

Insulin infusion mL/hr

Insulin concentration U/mL

Nutrition

☒ Enteral

Type

Rate mL/hr

☒ Variable ☐ Constant

☒ Parenteral

Type

Rate mL/hr

☐ Maintenance

Cancel OK

17:36

Patient Data is entered
into STAR, creating a
new current patient

PATIENT OVERVIEW

Patient Name
and Hospital ID

Last Blood Glucose
Measurement

Current Nutrition
Treatment

STAR_android

Patient
SMITH, John **BOND007**

Actions

Calculate treatment

Make changes **Stop STAR**

Details

BG **14.0 mmol/L at 17:38**
Next BG at: 18:38. BG target range: 4.4 - 8.0 mmol/L.

Insulin

Bolus **6.0mL at 17:38**

Infusion **0.0mL/hr**
Variable

Nutrition

Enteral **60 mL/hr**
Variable (91% goal)

Parenteral **20 mL/hr**
Constant

Maintenance **none**
Constant

Dextrose bolus **0.0mL**


Select patient **Patient history**

17:43

Current Insulin
Treatment

CALCULATING A TREATMENT

Another Blood Glucose measurement is to be taken at the prescribed time

STAR_android

Patient
SMITH, John **BOND007**

Actions

Calculate treatment

Make changes **Stop STAR**

Details

Blood Glucose **14.0 mmol/L at 17:38**
Next BG at: 18:38. BG target range: 4.4 - 8.0 mmol/L.

Insulin

Bolus **6.0mL at 17:38**
Infusion **0.0mL/hr**
Variable

Nutrition

Enteral **60 mL/hr**
Variable (91% goal)
Parenteral **20 mL/hr**
Constant
Maintenance **none**
Constant
Dextrose bolus **0.0mL**

Select patient Patient history

17:43

Having taken this Blood Glucose measurement press 'Calculate Treatment'

ENTERING A BLOOD GLUCOSE MEASUREMENT

Enter Blood Glucose measurement taken

The screenshot shows the STAR_android app interface. At the top, it says "Blood glucose" and "Please enter the new blood glucose measurement". Below this is a numeric keypad with buttons for digits 0-9 and a "Del" button. The value "10.0" is entered in the input field, followed by "mmol/L". Below the keypad is a date and time selection area with buttons for "Nov", "20", "2011", "18", and "38". At the bottom, there are "Back" and "OK" buttons. The status bar at the very bottom shows the time "17:44".


BG entered as '1' '0' '0' - Key pad automatically inserts decimal place.

Time displayed is current time, this can be edited

CALCULATED TREATMENT RECOMMENDATION

Historical BG
plot is in **RED**

Recommended
Insulin and
Nutrition
treatments for
the next hour

The **BLUE**
bands show the
90%
confidence
interval for the
patients future
Blood Glucose


UPDATED INFORMATION

The Patient
Overview Screen
has been updated

Last Blood
Glucose
Measurement

Current Nutrition
Treatment

STAR_android

Patient
SMITH, John **BOND007**

Actions
Calculate treatment
Make changes **Stop STAR**

Details
BG **10.0 mmol/L at 18:38**
Next BG at: 19:38. BG target range: 4.4 - 8.0 mmol/L.

Insulin
Bolus **6.0mL at 18:38**
Infusion **0.0mL/hr**
Variable


Nutrition
Enteral **43 mL/hr**
Variable (65% goal)
Parenteral **20 mL/hr**
Constant
Maintenance **none**
Constant
Dextrose bolus **0.0mL**

Select patient Patient history

Current Insulin
Treatment

STANDARD PROCEDURE

STAR_android

Patient
SMITH, John **BOND007**

Actions
Calculate treatment

Make changes Stop STAR

Details

BG 10.0 mmol/L at 18:38
Next BG at: 19:38, BG target range: 4.4 - 8.0 mmol/L

Insulin Bolus 6.0mL at 18:38
Infusion 0.0mL/hr Variable

Nutrition Enteral 43 mL/hr Variable (65% goal)
Parenteral 20 mL/hr Constant
Maintenance none Constant
Dextrose bolus 0.0mL

Select patient Patient history

The procedure is the same each time a Blood Glucose Measurement is taken

STAR_android

Blood glucose

Please enter the new blood glucose measurement

5.5 mmol/L

7 8 9
4 5 6
1 2 3
0 Del

+ + + + +
Nov 20 2011 22 38
- - - - -


Training mode is active: time set to the next scheduled BG measurement

Back OK

MULTIPLE TREATMENT OPTIONS


If the patients Blood Glucose is lower than 8.0 mmol/L longer treatment options may be offered

The best treatment option is selected, based on the users discretion and perception of patient stability

The **BLUE** bands show the 90% confidence interval for future Blood Glucose for each treatment option


HYPOGLYCAEMIC EVENT

In the event of a low Blood Glucose event, a warning screen will appear


HYPOGLYCAEMIC EVENT

Insulin is to be turned off,
and additional nutrition
recommended

RESULTS

Very good glycemic control
of a 'difficult' patient

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